OpaCement OpaCement

With our continuous research and experience, we have been able to develop the new upgraded PMMA Bone Cement OpaCem.

The OpaCem cement includes 60% of X-ray contrast agent Zirconium Dioxide which allows to visualize the cement under fluoroscopy more clearly and determine the exact placement of cement.





270, Osongsaengyeong 1-ro, Osong-eup, Heungdeok-gu, Cheonju-si, Chungcheongbuk-do, KOREA
Tel: +82-43-230-8914 | Fax: +82-43-217-1988

E-mail: madina@meta-biomed.com




Indications

Severely painful osteoporotic compression fractures Painful vertebral body tumors(metastasis or myeloma) Symptomatic Vertebral Hemangioma


Features

High radiopaque cement includes 60% of ZrO2 Optimized viscosity for Vertebroplasty and Kyphoplasty Short mixing and preparing time allows to inject the cement immediately Good mechanical properties

Chart of handling time at Various Temperature



High radiopacity property



Pic 1 shows the Al step-wedge by 3, 4, 5, 6mm and below Pic shows the brightness of each product which can be compared to Al step-wedge The following table shows the radiopacity property of OpaCem cement, NTCem Spine(old version), and two other brands.

	No	Product Name	Contrast agent name and amount added	AI(mm)
	, 1	NTCem Spine	ZrO2 30%	2.2
V	2	OpaCem	ZrO2 60%	4.6
-	3	"B" brand	ZrO2 60%	4.7
	4	"S" brand	BaSO4 30%	2.4



Product Specification

Odel Name	Order No	Components	Weight(g)
		Powder	20g
OpaCem	BG001	Liquid	8 . 5g